

CLEVER CHOICE™ mini

BLOOD GLUCOSE MONITORING SYSTEM





Quick Start User Guide

Dear CLEVER CHOICE mini System Owner:

Thank you for choosing CLEVER CHOICE mini Blood Glucose Monitoring System to take control of your diabetes.

CLEVER CHOICE mini will help you keep track of your blood glucose levels.

The Quick Start User Guide will assist you in setting up the system and testing your blood glucose levels. Please read it carefully before use. For complete information on the system's use and features, please refer to the Owner's Manual.



1. Set The Meter

1.1 Press SET to start setting. Year flashes first.

1.2 Press **M** to change values. Hold down **M** to increase the number faster. Press SET to move on to the next step.

1.3 Date flashes. Follow **1.2** to make changes. The sequence of settings is:
 (1) the date
 (2) 12h/ 24h
 (3) time
 (4) memory deletion
 (5) alarm function → OFF

The meter will automatically turn off when the setting is complete. For memory deletion, please refer to the Owner's Manual for details.

2. Set The Lancing Device

2.1 Remove the cap of the lancing device.

2.2 Firmly insert a lancet into the lancet holder until fully in place.

2.3 Remove the protective disk.

2.4 Replace the lancing device cap. For AST, please replace the clear cap instead.

2.5 Select the depth of penetration you require: 1-2 for soft or thin skin, 3 for average skin, 4-5 for thick or calloused skin.

2.6 Slide the cocking control back until it clicks. Set aside for later use.

3. Test Your Blood Glucose

3.1 Tear off the security label (Tamper Seal). When opening the vial for the first time, record the date on the label.

3.2 Wash and dry your hands thoroughly. Take a test strip out and then close the vial tightly.

3.3 Insert a test strip to turn on the meter. The sequences of display are: "CHK" and " " → flashing " " with date and time.

3.4 Take the pre-set lancing device. Hold the lancing device firmly against the puncture site. Press the release button to obtain a drop of blood.

3. Test Your Blood Glucose

3.5 Gently squeeze around the punctured area to obtain enough blood.

3.6 While " " is flashing on the display, hold and touch the absorbent hole of the test strip to the drop of blood until the confirmation window is completely filled.

3.7 After the meter counts to 0, your test result will appear along with the date and time and will be automatically stored in the memory.

3.8 Eject the used test strip. Dispose of the used strip and lancet carefully.

For testing with control solution, please refer to the Owner's Manual for detailed information.

4. Read Your Result

Your blood glucose result will display the following messages:

Lo	Hi
Appears when glucose < 20 mg/dL (1.1mmol/L)	Appears when glucose > 600 mg/dL (33.3mmol/L)

Reference values:

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meals	Less than 100 mg/dL (5.6 mmol/L)
2 hours after a meal	Less than 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2008). Clinical Practice Recommendations. Diabetes Care, 31 (Supplement 1): S1-108.

Please consult your doctor to determine a target range that works best for you.

Distributed by Simple Diagnostics P.O. Box 128, Williston Park, NY 11596 USA 24-HOUR TOLL-FREE TELEPHONE ASSISTANCE 1-877-DIABETK 1-877-342-2385 www.SimpleDiagnostics.com

For in vitro diagnostic use only. Read instructions before use. For self-testing. 315-42610-001

If your lancing device differs from the one shown here, please refer to the manufacturer's manual to ensure proper usage.