Limited Two-Year Warranty

Simple Diagnostics guarantees that its Digital Blood Pressure Monitor will be free from manufacturing defects under normal use for as long as the original purchaser owns the product. This warranty covers only normal use, and does not apply to use in any clinical or commercial applications. This warranty does not cover batteries or other power sources that may be provided with or used with the Digital Blood Pressure Monitor. This warranty is voided if the Digital Blood Pressure Monitor is misused or abused in any manner.

If the Digital Blood Pressure Monitor fails to operate during the time the original purchaser owns it, return it postage prepaid along with \$5.00 for return shipping and insurance to Simple Diagnostics, Attn: Repair Department, PO Box 128, Williston Park, NY 11596. When returning a product, please also include a letter with your name, address, phone number and brief description of the specific problem. Please carefully package the product to avoid any damage that may occur while in the transit. Packaging insurance with returned receipt is recommended. Simple will repair or replace the defective unit at the option of Simple Diagnostics. Repair or replacement of the defective unit, at the warrantor's option, is the sole remedy under this warranty.

ANY IMPLIED WARRANTIES WHICH THE PURCHASER MAY HAVE ARE LIMITED IN DURATION TO THE TIME THAT THE ORIGINAL CONSUMER PURCHASER OWNS THE PRODUCT. Some states do not allow limitations or how long an implied warranty lasts, so the above limitation may not apply to you. This warranty constitutes the warrantor's only responsibility and obligation to repair and/or replace materials or components, or refund the purchase price. **Simple will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others.** Simple shall in no event be liable to the purchaser for any amount in excess of the cost of repair and/or replacement of the unit, or the purchase price of the unit. Some states do not allow the exclusive or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

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Simple Diagnostics

Blood pressure is measured in millimeters (mm) of mercury (Hg) and is generally recorded with the systolic pressure (120) listed first and the diastolic pressure (80) listed second. The numbers are typically separated by a slash mark (/) as shown.

Both pressure readings, the SYSTOLIC and DIASTOLIC, are necessary for a physician to evaluate the status of a patient's blood pressure.

Please contact your physician for specific information regarding your own blood pressure.

What Influences Blood Pressure?

Many factors such as genetics, age, sex, altitude, physical activity, anxiety, muscular development, certain medications or even the time of day can influence blood pressure. Influences such as sleep or relaxation decrease blood pressure, while anxiety or exercise increase blood pressure.

WHO Blood Pressure Classifications

Standards for assessment of high or low blood pressure have been established by the World Health Organization (WHO) as shown on the following chart:

CATEGORY	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	COLOR INDICATOR
Optimal	<120	<80	GREEN
Normal	120-129	80-84	GREEN
High - Normal	130-139	85-89	GREEN
Stage 1 - Mild	140-159	90-99	YELLOW
Borderline Hypertension	160-179	100-109	ORANGE
Stage 2 - Moderate Hypertension	≥180	≥110	RED
Stage 3 - Severe Hypertension	≥140	<90	YELLOW
Isolated Systolic Hypertension	140-159	<90	ORANGE
	≥160	<90	RED

This chart is only a general guideline. Contact your physician or trained healthcare professional to determine your NORMAL blood pressure.

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Index

Introduction	3
Product Identification & Caution	4
Why Monitor Your Blood Pressure at Home?	5
General Blood Pressure Information	5-7
Important Information Before Use	7
World Health Organization (WHO) Indicator	8
Irregular Heartbeat Detection	9
Battery Installation/Replacement	10
Setting the Date/Time	11
Voice Language Setting	12
Applying Your Blood Pressure Cuff	13
Taking Your Blood Pressure Reading	14-15
Recalling Measurements in Memory	15
Description of Display Symbols	16
Troubleshooting	16
Care and Maintenance	17
Product Specifications	17
Sample Blood Pressure Log	18

Introduction

To achieve the maximum benefit from your blood pressure monitor, we recommend that you first consult with your physician or trained healthcare professional.

Thank you for purchasing a Simple Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

The method of measurement that your Automatic Monitor uses is called the oscillometric method. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor easy to use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Battery Installation/Replacement

1. Press down on tab and lift the cover.
 2. Insert or replace 4 "AA" batteries into the compartment, matching the indicated polarity symbols. Fig. 1.
 3. To close, match up and connect the bottom hooks of the cover, then push the top end.
- NOTE:** Replace the batteries when the Low Battery symbol appears on the display or when the display remains blank when the POWER button is pressed.
4. It is recommended to remove the batteries if the unit will not be used for an extended period of time.


Fig. 1

AC Power Adapter (optional)

An AC adapter is available separately for battery-free operation. The AC adapter plugs into the side of the unit as shown. Fig. 2. Please remember to practice electrical safety precautions while using the adapter.

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Fig. 2

Variations in Blood Pressure

Blood pressure is influenced by many factors and can change from moment to moment. Normally, blood pressure is lowest during sleeping periods and rises during the day. The graph below represents variations in blood pressure shown over a day with measurements taken every 5 minutes.

Important Information Before Use

1. Blood pressure readings should be interpreted by a physician or trained healthcare professional, who is familiar with your medical history.
2. Perform your measurement in a quiet place. You should be seated and relaxed.
3. Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
4. Rest 15 minutes prior to taking a reading.
5. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
6. Keep the monitor stable during measurements. Remain still; do not talk during measurements.
7. Record your daily readings on the chart provided.
8. Take your readings at the same time, each day (or as recommended by your physician).
9. Wait a minimum of 15 minutes between readings. The wait time may vary depending on individual physiological characteristics.
10. This device is intended for adult use.
11. The inflation or deflation process can be stopped by pressing the POWER button.

Taking Your Blood Pressure Reading

Proceed only after reading the previous sections of this manual.

1. Position the monitor on a flat, stable surface with the digital display panel in view. Insert the cuff tubing connector into the left side of the monitor. Fig. 1.
2. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart. Fig. 2. Relax your left hand.
3. Press the POWER button to turn the unit on and select the memory bank you want your reading to be stored in by pressing the MEMORY button to alternate between bank 1 or 2. Press the POWER button to confirm your selection.


Fig. 2

NOTA: The monitor will begin to guide you with voice prompts if the voice function is turned on. This unit will run a self-test, Fig. 3, and show the last measured reading, Fig. 4, or if there is no previous measurement, "0" as shown in Fig. 5.

5. When the display appears as shown, Fig. 6 the monitor is ready and will begin to guide you with voice prompts if the voice function is on. The cuff will automatically inflate to begin.
6. The display will show the increasing pressure in the cuff. Along with the measurement progress on the left side of the screen.
7. The cuff will automatically inflate to approximately 190 mmHg or 40 mmHg higher than the last systolic pressure stored in memory.
8. Once the cuff inflates to the appropriate pressure level, the measurement will begin.
9. The unit will begin to automatically deflate and the numbers appearing on the display will decrease.
10. When all the air is deflated, the measurement is complete. The unit will display your blood pressure (systolic and diastolic) measurements and pulse rate. Fig. 7.

With the voice function on (L1 or L2), the monitor will say the measurement with the language you selected according to the classification standard of WHO (World Health Organization).

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Setting the Date/Time

NOTE: When the unit is off, the date and time will still appear at the top of the screen.

1. When the unit is off, press and release both the POWER and MEMORY button at the same time. The month will begin to flash. Fig. 1.
2. Press the MEMORY button to change until the desired number appears.


Fig. 1

NOTE: To scroll through the numbers quicker, hold down the MEMORY button.

3. Press the POWER button to set the entry.

4. The date, Fig. 2, followed by the hour and minutes will blink. Repeat steps 2 and 3 to adjust the date/time setting.

5. After setting the minutes, the unit will automatically exit out of the date/time setting mode. Fig. 3

NOTE: The time is expressed in 24hr military time. The chart below explains the time conversion.

Regular Time (12 hours)	Military Time (24 hours)	Regular Time (12 hours)	Military Time (24 hours)
12 am	0	12 pm	12
1 am	1	1 pm	13
2 am	2	2 pm	14
3 am	3	3 pm	15
4 am	4	4 pm	16
5 am	5	5 pm	17
6 am	6	6 pm	18
7 am	7	7 pm	19
8 am	8	8 pm	20
9 am	9	9 pm	21
10 am	10	10 pm	22
11 am	11	11 pm	23

Fig. 3

Sample Blood Pressure Log


Name: _____ Age: _____ Weight: _____

Date	1/2	1/2	1/3
Time	7:00 AM	9:30 AM	7:30 AM
mmHg	130	125	120
	84	90	78
Pulse	68	71	69
Body Condition	Headache in AM		Tired, Restless Night

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INSTRUCTION MANUAL

Talking Monitor
 Premium Digital Blood Pressure Arm Monitor


SDI-886AT

Please read this instruction manual completely before operating this unit.

DETAILED GUIDEBOOK
 English • Spanish

Product Identification & Caution

CAUTION: You can stop the inflation or deflation process anytime by pressing the Power button.

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Why Monitor Your Blood Pressure at Home?

Many people experience "White Coat Syndrome", which is a temporary elevation of blood pressure during a visit at their physician's office. Although only temporary, it can cause blood pressure to rise and give a false indication of a person's true blood pressure.

Whether or not you experience "White Coat Syndrome", home blood pressure monitoring provides:

- The opportunity to supplement your physician's office measurements.
- An accurate indication of any changes.
- A record to assist your physician in evaluating your health and making important decisions.
- Consistent, daily measurements of your blood pressure.

NOTE: Variations in your individual readings should only be interpreted by your physician or trained healthcare professional.

General Blood Pressure Information

What is Blood Pressure?
Blood pressure is the pressure that is exerted by blood flowing against the walls of the blood vessels throughout your body.

Your heart, which is the center of the circulatory system, provides the force for the blood to flow or circulate. When your heart contracts or beats, the blood is forced through the blood vessels increasing the pressure. This is the highest pressure in the cycle or what is referred to as SYSTOLIC blood pressure. In between beats, your heart relaxes and your blood pressure decreases. This is referred to as DIASTOLIC blood pressure.

This complete series of events, which occurs in a single heartbeat, is known as the CARDIAC CYCLE.

Your Simple Automatic Blood Pressure Monitor will automatically read your blood pressure and display both systolic and diastolic readings on the screen. Your systolic will be positioned as the upper number and the diastolic reading will be the lower number.

(systolic)
120/80
(diastolic)

World Health Organization (WHO) Indicator

The WHO Indicator, located on the left side of the display, will identify your reading within a range in accordance with the guidelines set by the World Health Organization (WHO). See the chart below for reference.

SYS	DIA	Stage 3 Hypertension
SYS	DIA	Systolic pressure: ≥ 180 mmHg Diastolic pressure: ≥ 110 mmHg
SYS	DIA	Systolic pressure: 160-179 mmHg Diastolic pressure: 100-109 mmHg
SYS	DIA	Systolic pressure: 140-159 mmHg Diastolic pressure: 90-99 mmHg
SYS	DIA	High Normal Systolic pressure: 130-139 mmHg Diastolic pressure: 85-89 mmHg
SYS	DIA	Normal Systolic pressure: 120-129 mmHg Diastolic pressure: 80-84 mmHg
SYS	DIA	Optimal Systolic pressure: < 120 mmHg Diastolic pressure: < 80 mmHg

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Irregular Heartbeat Detection

Your digital blood pressure monitor features Irregular Heartbeat Detection. This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur.

When an irregular heartbeat is detected, the IHB icon will appear on the display. Fig. 1.

NOTE: Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently at the end of a reading.

Voice Language Setting

To set the desired voice language:

- When the unit is off, you can select the desired voice language by pressing and holding the MEMORY button.
- Now the screen will cycle through L0, L1 and L2 as shown in Fig. 1.
L0 - voice function is off
L1 - voice is in English
L2 - voice is in Spanish
- The monitor will say the corresponding language with the displayed language code.
- Once the desired language has been selected, release the MEMORY button.
- After the desired language has been selected, you can adjust the volume by pressing and holding down the MEMORY button until the bars on the right of the screen begin to increase Fig. 2.
- Continue to hold the MEMORY button down until the volume reaches the desired level, then release the MEMORY button to set.
- The unit will automatically turn off after the desired volume level has been set.

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Applying Your Blood Pressure Cuff

Before applying your blood pressure cuff, be sure you have selected the appropriate cuff.
Small Adult Cuff fits arm circumference 8-5/8" - 11-3/4" (22-30 cm)
Adult Cuff fits arm circumference 11-3/4" - 16-1/2" (30-42 cm)
Large Adult Cuff fits arm circumference 16-1/2" - 18-7/8" (42-48 cm)

NOTE: If for any reason you are unable to or should not use your left arm, please modify the following instructions and apply the cuff to your right arm.

- Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
- Be seated with your feet flat on the floor.
- Position the cuff on a solid surface with the tubing facing UP and AWAY from you. Fig. 1. The metal bar on the cuff should be to the left of the tubing.
- Widen or open the cuff by pulling or rolling the bottom of the cuff towards the right. Fig. 2. This should open the cuff, creating a cylinder. Do not extend the cuff beyond the metal bar.
- Insert your arm into the cuff (cylinder). Position the (Q) mark over the main artery (on the inside of your arm). Fig. 3.
- The bottom edge of the cuff should be positioned approximately one inch ABOVE the elbow joint. Fig. 4.
- Reaching underneath your left arm with your right hand, pull the end of the cuff towards your body to tighten the cuff. Fig. 5. Wrap and secure the cuff making sure that the (Q) mark remains as shown. Fig. 3.
- The cuff should fit comfortably, yet snugly around your arm. You should be able to insert two fingers between your arm and the cuff.

NOTE: The cuff size is suitable for use when the vertical "Index" mark is within the horizontal "OK" range. A different sized cuff is needed if the "Index" mark is outside of the "OK" range.

Description of Display Symbols

SYS	Indicates a measurement within Stage 3 Hypertension.	Shows the pulse rate per minute.
DIA	Indicates a measurement within Stage 2 Hypertension.	Occurs when an error was made during measurement.
SYS	Indicates a measurement within Stage 2 Hypertension.	Appears when there is an irregular heartbeat detected during measurement.
DIA	Indicates a measurement within Stage 1 Hypertension.	Appears when the cuff is unstable or it has too much air in it. Once the arrow stops flashing, the unit is ready to take a measurement.
SYS	Indicates a measurement within Stage 1 Hypertension.	Appears when a measurement error has occurred. Restart the measurement.
DIA	Appears when batteries should be replaced.	Appears when recalling a measurement from memory.
		Indicates the memory user's memory bank.

Troubleshooting

If any abnormality occurs during use, please check and correct the following:

Condition	Correction
Display is blank when power is on.	Check and correct the polarity of the installed batteries. Reinstall or replace batteries.
Measurement incomplete or abnormally low or high values displayed.	Review and follow "Applying Your Blood Pressure Cuff" and "Taking Your Blood Pressure Reading" sections. Replace the batteries.
EE symbol displays.	Do not talk or move during the measurement. Review "Applying Your Blood Pressure Cuff" section.
Measurements are different from those typically measured by physician or every measurement reading is different.	Remember, blood pressure readings are influenced by physical and mental conditions and/or even the time of day. Daily measurements should be interpreted by your physician.

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CLEVER CHOICE™ PO Box 128, Williston Park, NY 11596

Care and Maintenance

- Disconnect the cuff tubing from the monitor and remove batteries prior to storing.
- Only use a soft, dry cloth to clean your blood pressure monitor.
- Avoid using any types of liquids on the monitor or cuff.
- Do not store the unit where it will be exposed to direct sunlight, dust or humidity.
- Avoid extreme temperatures.
- Never disassemble the monitor or cuff.
- Dropping or subjecting your blood pressure monitor to strong shocks should be avoided.

Product Specifications

Name and model number	Talking Monitor-Premium Digital Blood Pressure Arm Monitor, SDI-886AT
Display system	Digital display/LCD
Measuring method	Oscillometric
Power source	4 "AA" batteries
Measuring range	Pressure: 0-300 mmHg Pulse: 30-180 beats/minute
Accuracy	Pressure: ± 3 mmHg Pulse: ± 5 of reading
Inflation	Pump Drive
Deflation	Automatic Pressure Release Valve
Memory	Built-in memory enabling display of up to 120 (60 per memory bank) measurements
Automatic shut-off	Approximately 1 minute after last button operation
Operation environment	Temperature: 41°F - 104°F (5°C - 40°C) Humidity: $< 85\%$
Storage environment	Temperature: 4°F - 131°F (20°C - 55°C) Humidity: $< 95\%$
Monitor dimensions	4-1/4" x 6" x 2-3/8"
Weight	10.5 oz. (Without batteries)
Arm circumference range	Adult size cuff fits arm circumference: 11-3/4" - 16-1/2" (30-42 cm), Small adult size cuff fits arm circumference: 8-5/8" - 11-3/4" (22-30 cm), AC Adapter
Accessories	Cuff, detailed guidebook, quick start guide, blood pressure log, storage case
Options (sold separately)	Large adult size cuff fits arm circumference: 16-1/2" - 18-7/8" (42-48 cm), Small adult size cuff fits arm circumference: 8-5/8" - 11-3/4" (22-30 cm), AC Adapter

Specifications are subject to change without notice